



Product Datasheet

UBQLN4/CIP75/Ubiquilin 4 Antibody (orb1519632)

Description	UBQLN4/CIP75/Ubiquilin 4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	residues 2-601 (full length)
Target	UBQLN4/CIP75/Ubiquilin 4
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	50 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	Q9NRR5
NCBI	NP_064516.2
Dilution Range	WB - 1:1,000; IP - 20 µl/1 mg lysate; IHC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.; ICC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections.
Expiration Date	12 months from date of receipt.



Detection of human UBQLN4 by western blotting



Detection of human UBQLN4/CIP75/Ubiquilin 4 by flow cytometry