

Product Datasheet

HSP22 antibody [Out of stock] (orb151974)

Description

Rabbit polyclonal to Hsp22 (ATTO-594).

protein that is a member of...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ICC, IHC, IP, WB

Applications
Immunogen

Human HSP22

Target

HSP22

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-181 was sufficient for detection of HSP22 in 10 µg of rat tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

22kDa

Uniprot ID
[Q9UJY1](#)
NCBI
[NP_055180.1](#)
Entrez
[26353](#)
Dilution Range

WB (1:1000), ICC/IF (1:100), IHC (1:100)

Expiration Date

12 months from date of receipt.