

Product Datasheet

MEK1 Antibody (orb1519758)

Description	MEK1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	between 1 and 50
Target	MEK1
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	Q02750
NCBI	NP_002746.1
Dilution Range	WB - 1:1000; IP - 20 µl/1 mg lysate
Expiration Date	12 months from date of receipt.