

Product Datasheet

HSP22 antibody (orb151981)

Description	Rabbit polyclonal to Hsp22 (Biotin). Hsp22 (HSPB8) is a 196-amino acid protein that is a member of t...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC, IP, WB
Immunogen	Human HSP22
Target	HSP22
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-181 was sufficient for detection of HSP22 in 10 µg of rat tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	22kDa
Uniprot ID	Q9UJY1
NCBI	NP_055180.1
Entrez	26353
Dilution Range	WB (1:1000), ICC/IF (1:100), IHC (1:100)
Expiration Date	12 months from date of receipt.