

Product Datasheet

Calnexin-CT antibody [Out of stock] (orb151989)

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Reactivity	Aves, Bovine, Canine, Drosophila, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	ATTO-488
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	Dog Calnexin C-terminal synthetic peptide conjugated to KLH. Identical to human, mouse and rat calnexin sequences over these residues.
Target	Calnexin
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:2000 dilution of SPC-182 was sufficient for detection of Calnexin in 10 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	90kDa
Uniprot ID	P24643
NCBI	NP_001003232.1
Entrez	403908
Dilution Range	WB (1:2000), ICC/IF (1:100), IHC (1:100)
Expiration Date	12 months from date of receipt.