

Product Datasheet

Calnexin-CT antibody (orb151998)

Description

Rabbit polyclonal antibody against Calnexin

Species/Host	Rabbit
Reactivity	Aves, Bovine, Canine, Drosophila, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Biotin
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	Dog Calnexin C-terminal synthetic peptide conjugated to KLH. Identical to human, mouse and rat calnexin sequences over these residues.
Target	Calnexin
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:2000 dilution of SPC-182 was sufficient for detection of Calnexin in 10 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	90kDa
Uniprot ID	P24643
NCBI	NP_001003232.1
Entrez	403908
Dilution Range	WB (1:2000), ICC/IF (1:100), IHC (1:100)
Expiration Date	12 months from date of receipt.