

Product Datasheet

BRD2 Antibody (orb1519984)

Description	BRD2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	between 675 and 725
Target	BRD2
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	40 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P25440
NCBI	NP_005095.1
Dilution Range	WB - 1:1000; IP - 20 µl/1 mg lysate
Expiration Date	12 months from date of receipt.