


Product Datasheet

Cytokeratin Antibody (orb1520125)


Description	Cytokeratin Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Human epidermal callus
Target	Cytokeratin
Preservatives	Phosphate Buffered Saline (PBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Monoclonal
Uniprot ID	P08779, P19012, P08727, P13645, P02533
NCBI	NP_002266, NP_005548, NP_000412, NP_002267, NP_000517
Dilution Range	IHC - 1:100 - 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.; IHC-IF - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.; ICC-IF - 1:100 - 1:500
Expiration Date	12 months from date of receipt.



Detection of human Cytokeratin by immuno...



Detection of human Cytokeratin by immuno...



Detection of human Cytokeratin by immuno...