

Product Datasheet

ACE2(18-740)/Rb-Fc Antibody (orb1520211)

Description	ACE2(18-740)/Rb-Fc Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	HEK293-derived human ACE2 Gln18-Ser740, with a C-terminal Myc/TEV/HIS6/Rb-Fc Tag
Target	ACE2(18-740)/Rb-Fc
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BYF1
NCBI	NP_068576.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µl/mg lysate
Expiration Date	12 months from date of receipt.