



Product Datasheet

FAT4 Antibody (orb1520253)

Description	FAT4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 4931 and 4981
Target	FAT4
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6V0I7
NCBI	Q6V0I7.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 8-15 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human FAT4 by western blot....



Detection of human FAT4 by western blot ...