



Product Datasheet

TTC37 Antibody (orb1520304)

Description	TTC37 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 1 and 50
Target	TTC37
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6PGP7
NCBI	Q6PGP7.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 20-40 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of
human and
mouse TTC37
by we...



Detection of
human TTC37
by western
blot...