



Product Datasheet

METTL14 Antibody (orb1520323)


Description	METTL14 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 25 and 75
Target	METTL14
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9HCE5
NCBI	Q9HCE5.2
Dilution Range	WB - 1:10,000 - 1:25,000; IP - 50-100 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human METTL14 by western bl...



Detection of mouse METTL14 by western bl...



Detection of human METTL14 by western bl...