


Product Datasheet

PPP1R21 Antibody (orb1520367)

Description	PPP1R21 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 730 and 780
Target	PPP1R21
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6ZMI0
NCBI	Q6ZMI0.1
Dilution Range	IP - 15-20 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human PPP1R21 by western bl...