



Product Datasheet

HEATR3 Antibody (orb1520400)

Description	HEATR3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 25 and 75
Target	HEATR3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7Z4Q2
NCBI	Q7Z4Q2.2
Dilution Range	WB - 1:10,000 - 1:25,000; IP - 50-100 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse HEATR3 by W...



Detection of human HEATR3 by western blo...