


Product Datasheet

FAM134A Antibody (orb1520416)

Description	FAM134A Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 325 and 375
Target	FAM134A
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8NC44
NCBI	Q8NC44.3
Dilution Range	IP - 15-20 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human FAM134A by western blotting.