



Product Datasheet

alpha-Tubulin Antibody (orb1520421)


Description	alpha-Tubulin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 175 and 225
Target	alpha-Tubulin
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q71U36
NCBI	NP_001257328.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2-10 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse TUBA1A by western blo...



Detection of human TUBA1A by western blo...



Detection of human TUBA1A by western blo...