


Product Datasheet

LPCAT4 Antibody (orb1520458)

Description	LPCAT4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 474 and 524
Target	LPCAT4
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q643R3
NCBI	Q643R3.1
Dilution Range	IP - 8 - 15 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human LPCAT4 by western blotting