


Product Datasheet


MAPK1IP1L Antibody (orb1520517)

Description	MAPK1IP1L Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 195 and 245
Target	MAPK1IP1L
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8NDC0
NCBI	Q8NDC0.4
Dilution Range	WB - 1:1,000 - 1:5,000; IP - 3 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MAPK1IP1L by western

...



Detection of human MAPK1IP1L by western

...