


Product Datasheet

C19ORF33 Antibody (orb1520526)

Description	C19ORF33 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 1 and 50
Target	C19ORF33
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9GZP8
NCBI	Q9GZP8.1
Dilution Range	WB - 1:2,000 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of human and mouse C19ORF33 by...