


Product Datasheet

TTYH3 Antibody (orb1520585)

Description	TTYH3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 425 and 475
Target	TTYH3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9C0H2
NCBI	Q9C0H2.3
Dilution Range	WB - 1:10,000 - 1:25,000
Expiration Date	12 months from date of receipt.



Detection of human TTYH3 by western blot...