


Product Datasheet

PDE12 Antibody (orb1520709)

Description	PDE12 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 559 and 609
Target	PDE12
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6L8Q7
NCBI	Q6L8Q7.2
Dilution Range	WB - 1:2,000 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of human PDE12 by western blot...