



Product Datasheet

SF3B5 Antibody (orb1520805)

Description	SF3B5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1 and 50
Target	SF3B5
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BWJ5
NCBI	Q9BWJ5.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 20 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse SF3B5 by we...



Detection of human SF3B5 by western blot...