



Product Datasheet

METAP1 Antibody (orb1520808)


Description	METAP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 336 and 386
Target	METAP1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P53582
NCBI	P53582.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 20 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse METAP1 by western blotting



Detection of human METAP1 by western blotting



Detection of human METAP1 by western blotting