

Product Datasheet

**HSP90 (P. falciparum) antibody [Out of stock]
(orb152081)**

Description

Rabbit polyclonal to Hsp90 antibody

Species/Host	Rabbit
Conjugation	AP
Tested Applications	ICC, IP, WB
Immunogen	Recombinant full length PfHSP90
Target	HSP90
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1.56 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.7 µg/ml of SPC-187 was sufficient for detection of PfHSP90 in 20 µg of P. falciparum lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	86kDa
Uniprot ID	Q8IL32
NCBI	XP_001348591.1
Entrez	811999
Dilution Range	WB (1:2000), ICC/IF (1:50)
Expiration Date	12 months from date of receipt.