

Product Datasheet

Cpn10 antibody [Out of stock] (orb152093)

Description

Rabbit polyclonal to Cpn10 (ATTO-594).

Species/Host

Rabbit

Reactivity

Bovine, Canine, Frog, Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Sheep

Conjugation

ATTO-594

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Human Cpn10 peptide AA 91-101

Target

Cpn10

Preservatives

PBS pH 7.3, 0.02% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-193 was sufficient for detection of Cpn10 in 10 µg of Hela lysates by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

10kDa

Uniprot ID
P61604
NCBI
NP_002148.1
Entrez
3336
Dilution Range

WB (1:1000), IHC (1:100)

Expiration Date

12 months from date of receipt.