

Product Datasheet

HSP110 antibody (orb152125)

Description

Rabbit polyclonal to Hsp110 (ATTO-488).

	proteins known as...
Species/Host	Rabbit
Reactivity	Bovine, Fish, Hamster, Human, Monkey, Mouse, Rat, Sheep, Yeast
Conjugation	ATTO-488
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide derived from the sequence of hamster HSP110; sequence identical to human and mouse
Target	HSP110
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-195 was sufficient for detection of HSP110 in 10 µg of human cell line mixed lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	110kDa
Uniprot ID	Q92599
NCBI	NP_006635.3
Entrez	10808
Dilution Range	WB (1:1000), IHC (1:100)
Expiration Date	12 months from date of receipt.