

Product Datasheet

PDIA4/ERP72 Antibody (orb1521408)

Description	PDIA4/ERP72 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 595 and 645
Target	PDIA4/ERP72
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P13667
NCBI	P13667.2
Dilution Range	WB - 1:10,000 - 1:25,000; IP - Not recommended
Expiration Date	12 months from date of receipt.