

Product Datasheet

ITGA2/CD49B/Integrin alpha-2 Antibody (orb1521516)

Description	ITGA2/CD49B/Integrin alpha-2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 25 and 75
Target	ITGA2/CD49B/Integrin alpha-2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P17301
NCBI	P17301.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended
Expiration Date	12 months from date of receipt.