

Product Datasheet

Piccolo antibody (orb152168)

Description	Rabbit polyclonal to Piccolo (Biotin). Piccolo, also referred to as Aczonin, is a large protein whic...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	Recombinant Full length Protein, BC001304.1
Target	Piccolo
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-197 was sufficient for detection of Piccolo in 10 µg of rat brain tissue lysate by colorimetric immunoblot analysis using goat ant rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	550kDa
Uniprot ID	Q9Y6V0
NCBI	NP_055325.2
Entrez	27445
Dilution Range	WB (1:1000), IHC (1:400)
Expiration Date	12 months from date of receipt.