

Product Datasheet

ATXN1L Antibody (orb1521681)

Description	ATXN1L Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 250 and 300
Target	ATXN1L
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P0C7T5
NCBI	P0C7T5.1
Dilution Range	WB - Not recommended. orb1521681 has not performed satisfactorily when used for WB of ATXN1L in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of ATXN1L.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.