

## Product Datasheet

### DIEXF/C1orf107 Antibody (orb1521688)

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | DIEXF/C1orf107 Antibody                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                    |
| <b>Reactivity</b>          | Human                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | IP                                                                                        |
| <b>Immunogen</b>           | Between 150 and 200                                                                       |
| <b>Target</b>              | DIEXF/C1orf107                                                                            |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide                    |
| <b>Form/Appearance</b>     | Whole IgG                                                                                 |
| <b>Concentration</b>       | 1000 µg/ml                                                                                |
| <b>Storage</b>             | 2 - 8°C                                                                                   |
| <b>Note</b>                | For research use only                                                                     |
| <b>Application notes</b>   | Format: Whole IgG                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q68CQ4</a>                                                                    |
| <b>NCBI</b>                | <a href="#">Q68CQ4.2</a>                                                                  |
| <b>Dilution Range</b>      | WB - Not recommended. Use rabbit anti-DIEXF antibody A305-122A.; IP - 2 - 10 µg/mg lysate |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |