

Product Datasheet

Bassoon antibody [Out of stock] (orb152176)

Description

Rabbit polyclonal to Bassoon (ATTO-488).

localized at the pre...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-488

Tested

IHC, WB

Applications
Immunogen

NM_003458.3 (AA 786-1041) N-terminal his-tagged fusion protein

Target

Bassoon

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-198 was sufficient for detection of Bassoon in 10 µg of rat brain tissue lysate by colorimetric immunoblot analysis using goat ant rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

420kDa

Uniprot ID
[Q9UPA5](#)
NCBI
[NP_003449.2](#)
Entrez
[8927](#)
Dilution Range

WB (1:1000), IHC (1:400)

Expiration Date

12 months from date of receipt.