

Product Datasheet

Bassoon antibody (orb152185)

Description	Rabbit polyclonal to Bassoon (Biotin). Bassoon is a 420 kDa protein that is localized at the presy...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	NM_003458.3 (AA 786-1041) N-terminal his-tagged fusion protein
Target	Bassoon
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-198 was sufficient for detection of Bassoon in 10 µg of rat brain tissue lysate by colorimetric immunoblot analysis using goat ant rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	420kDa
Uniprot ID	Q9UPA5
NCBI	NP_003449.2
Entrez	8927
Dilution Range	WB (1:1000), IHC (1:400)
Expiration Date	12 months from date of receipt.