

Product Datasheet

TLR4 antibody [Out of stock] (orb152192)

Description

Rabbit polyclonal to TLR4 (ATTO-390). The

comprises a fami...

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

ATTO-390

Tested

FACS, FC, ICC, IHC, WB

Applications
Immunogen

Developed against a synthetic peptide corresponding to amino acids 420-435 of human TLR4

Target

TLR4

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

2 µg/ml of SPC-200 was sufficient for detection of TLR4 in 100 ng of partial recombinant mouse TLR4 protein by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

75-80kDa

Uniprot ID

[O00206](#)

NCBI

[NP_612564.1](#)

Entrez

[7099](#)

Dilution Range

WB (1:500), IHC (1:50)

Expiration Date

12 months from date of receipt.