

Product Datasheet

SLX4IP/C20orf94 Antibody (orb1521931)

Description	SLX4IP/C20orf94 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 175 and 225
Target	SLX4IP/C20orf94
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q5VYV7
NCBI	Q5VYV7.1
Dilution Range	WB - Not recommended. Use rabbit anti-SLX4IP/C20orf94 antibody A304-994A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.