



Product Datasheet

KLHL36 antibody (orb1522)

Description	KLHL36 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human KLHL36/C16orf44 (51-150/616 aa)
Target	KLHL36
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70 kDa
Uniprot ID	Q8N4N3
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human lung tissue.