

Product Datasheet

TLR4 antibody (orb152203)

Description

Rabbit polyclonal to TLR4 (FITC). The Toll-like

family o...

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

FITC

Tested

FACS, FC, ICC, IHC, WB

Applications
Immunogen

Developed against a synthetic peptide corresponding to amino acids 420-435 of human TLR4

Target

TLR4

Preservatives

PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

2 µg/ml of SPC-200 was sufficient for detection of TLR4 in 100 ng of partial recombinant mouse TLR4 protein by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

75-80kDa

Uniprot ID
[O00206](#)
NCBI
[NP_612564.1](#)
Entrez
[7099](#)
Dilution Range

WB (1:500), IHC (1:50)

Expiration Date

12 months from date of receipt.