


Product Datasheet

HSD17B4/MFE-2/17-beta-HSD4 Antibody (orb1522309)

Description	HSD17B4/MFE-2/17-beta-HSD4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 600 and 650
Target	HSD17B4/MFE-2/17-beta-HSD4
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P51659
NCBI	NP_000405.1
Dilution Range	WB - Not recommended. orb1522309 has not performed satisfactorily when used for WB of HSD17B4/MFE-2/17-beta-HSD4 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of HSD17B4/MFE-2/17.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human HSD17B4/MFE-2 by west...