



Product Datasheet

gamma Adaptin/AP1G1 Antibody (orb1522360)


Description	gamma Adaptin/AP1G1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1 and 50
Target	gamma Adaptin/AP1G1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O43747
NCBI	NP_001119.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse gamma Adaptin/AP1G1 b...



Detection of human gamma Adaptin/AP1G1 b...



Detection of human gamma Adaptin/AP1G1 b...