



## Product Datasheet

### C2CD2L/TMEM24 Antibody (orb1522374)

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | C2CD2L/TMEM24 Antibody                                                 |
| <b>Species/Host</b>        | Rabbit                                                                 |
| <b>Reactivity</b>          | Human, Mouse                                                           |
| <b>Conjugation</b>         | Unconjugated                                                           |
| <b>Tested Applications</b> | IP, WB                                                                 |
| <b>Immunogen</b>           | Between 656 and 706                                                    |
| <b>Target</b>              | C2CD2L/TMEM24                                                          |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| <b>Form/Appearance</b>     | Whole IgG                                                              |
| <b>Concentration</b>       | 1000 µg/ml                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                |
| <b>Note</b>                | For research use only                                                  |
| <b>Application notes</b>   | Format: Whole IgG                                                      |
| <b>Clonality</b>           | Polyclonal                                                             |
| <b>Uniprot ID</b>          | <b>O14523</b>                                                          |
| <b>NCBI</b>                | <b>NP_001277403.1</b>                                                  |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                        |



Detection of human and mouse C2CD2L/TMEM24



Detection of human C2CD2L/TMEM24 by west...