



Product Datasheet

SOD1 Antibody (orb152243)

Description	Rabbit Anti-Human SOD1 (UbetaB) Polyclonal
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Within 50 AA of N-terminal region of SOD1, unfolded beta barrel region
Target	SOD1
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1µg/ml of SPC-205 was sufficient for detection of the unfolded beta barrel of SOD1 by colorimetric dot blot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
Uniprot ID	P00441
NCBI	CAG46542
Entrez	6647
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western Blot analysis of human, rat, mou...