



Product Datasheet

Beta-2-adaptin Antibody (orb1522465)


Description	Beta-2-adaptin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 575 and 625
Target	Beta-2-adaptin
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P63010
NCBI	NP_001273.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse Beta-2-adaptin by wes...



Detection of human Beta-2-adaptin by wes...



Detection of
human Beta-
2-adaptin by
wes...