


Product Datasheet

C12orf45 Antibody (orb1522490)

Description	C12orf45 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 75 and 125
Target	C12orf45
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8N5I9
NCBI	NP_689531.2
Dilution Range	WB - Not recommended. Use rabbit anti-C12orf45 antibody A304-707A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human C12orf45 by western b...