


Product Datasheet

DMXL1 Antibody (orb1522525)

Description	DMXL1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 125 and 175
Target	DMXL1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9Y485
NCBI	NP_005500.4
Dilution Range	WB - Not recommended. orb1522525 has not performed satisfactorily when used for WB of DMXL1 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of DMXL1.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human DMXL1 by western blot...