


Product Datasheet

OVCOV1 Antibody (orb1522646)

Description	OVCOV1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 315 and 365
Target	OVCOV1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NQQ7
NCBI	NP_057029.8
Dilution Range	WB - 1:1,000 - 1:5,000; IP - Not recommended
Expiration Date	12 months from date of receipt.



Detection of human and mouse OVCOV1 by W...