

Product Datasheet

SOD1 (EDI) antibody (orb152276)

Description	Rabbit polyclonal to SOD1 (PerCP). Superoxide dismutase (SOD) is an endogenously produced intracellu...
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	PerCP
Tested Applications	ELISA, IHC, WB
Immunogen	C-terminal region of SOD1, exposed dimer interface (EDI)
Target	SOD1
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1µg/ml of SPC-206 was sufficient for detection of the exposed dimer interface of SOD1 by colorimetric dot blot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
Uniprot ID	P00441
NCBI	CAG46542
Entrez	6647
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.