



Product Datasheet

DGCR14 Antibody (orb1522762)

Description	DGCR14 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 426 and 476
Target	DGCR14
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96DF8
NCBI	NP_073210.1
Dilution Range	WB - 1:500 - 1:2,500; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse DGCR14 by W...



Detection of human DGCR14 by western blo...