


Product Datasheet

GRASP65 Antibody (orb1522770)

Description	GRASP65 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 375 and 425
Target	GRASP65
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BQQ3
NCBI	NP_114105.1
Dilution Range	WB - Not recommended. Use rabbit anti-GRASP65 antibody A304-563A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human GRASP65 by western bl...