


Product Datasheet

RLIP76 Antibody (orb1522892)

Description	RLIP76 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 600 and 650
Target	RLIP76
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q15311
NCBI	NP_006779.1
Dilution Range	WB - Not recommended. Use rabbit anti-RLIP76 antibody A304-498A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human RLIP76 by western blo...