


Product Datasheet

Protocadherin-7 Antibody (orb1522936)

Description	Protocadherin-7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 1019 and 1069
Target	Protocadherin-7
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O60245
NCBI	NP_002580.2
Dilution Range	WB - 1:500 - 1:2,500; IP - Not recommended
Expiration Date	12 months from date of receipt.



Detection of human and mouse Protocadher...